



FMPTS-25

Food Microbiology Proficiency Testing Scheme

FMPTS

Food Microbiology Proficiency Testing Scheme **Preparation Instruction sheet**

Sample ID	Protocol	Page		
<u>General Information</u>		3		
PT-FM-101	FM-1V-1M(10g)	4-6		
PT-FM-100				
PT-FM-098				
PT-FM-097				
PT-FM-095				
PT-FM-090				
<hr/>				
PT-FM-096			FM-2V-1M(50g)	7-8
PT-FM-094				
PT-FM-092				
<hr/>				
PT-FM-093	FM-2V	9-10		
PT-FM-091				

FMPTS

Food Microbiology Proficiency Testing Scheme

Preparation Instruction sheet

GENERAL INFORMATION

Receipt and Storage

- On receipt of the test material, record the date and store refrigerated at 2-8^oC until ready to test.
- The test material should be analyzed in accordance with the deadlines shown on the contract/registration form.
- The test materials represent 'real' food samples, which may or may not contain the target organism(s), at a range of inoculum levels. Test materials may also contain background flora.
- Please consider all dilution factors when calculating the level of the target organisms in the original sample Where relevant guidance is given in the instructions.
- The test materials should be treated in the same way as routine samples and should therefore be tested using your routine laboratory equipment, reagents and methods.
- Where the instruction mentions a **resuscitate** step, the lyophilized test material should be left at room temperature for 45 minutes before proceeding with the test.
- Do not sub-divide samples prior to dilution unless stated in the instructions.
- Please treat the sample as contaminated material and use aseptic technique throughout the test procedure.
- Please ensure samples are mixed well at all stages before proceeding with testing.
- Reconstitute the test material according to the appropriate instructions for each sample type.

Precautions

- Microbiological test materials contain viable micro-organisms and are supplied on the understanding that the purchaser has suitably competent and qualified personnel to handle them safely. Test materials must only be opened in a laboratory by qualified personnel.
- Refer to the Safety Data Sheet for information on the safe handling and disposal of the test materials.

Reporting Results

- All results should be submitted using the NLPTS reporting form.

- If you need any help at all, please do not hesitate to contact our support team using the Details below.



Tel: +962 6 5506539

Email: pt-coordinator@naratechlab.com

FMPTS

Food Microbiology Proficiency Testing Scheme


Preparation Instruction sheet

Protocol FM1V1M	Sample code:			Description: 1 vial + 10g food matrix (milk powder)	
	PT-FM-101	PT-FM-100	PT-FM-098		
	PT-FM-097	PT-FM-095	PT-FM-090		

Step 1: prepare 90ml of your chosen diluent.



Step 2: add 10 ml from the broth prepared in step 1 to the vial after aseptically remove cap and rubber stopper.



FMPTS

Food Microbiology Proficiency Testing Scheme

Preparation Instruction sheet

Step 3: replace the vial stopper and shake gently to dissolve, allow to resuscitate before proceeding (45 min resuscitation).



Step 4: weigh 10g milk powder in sterile container



Step 5: Reconstitute the milk powder with the remaining amount of diluent prepared in step 1.




FMPTS

Food Microbiology Proficiency Testing Scheme

Preparation Instruction sheet

Step 6: after resuscitation add the vial content prepared in step 3 to the reconstituted Matrix from step 5. back wash 2 or 3 times to ensure all the freeze-dried material is recovered from the vial. (This represents the first dilution)




Step 7: proceed with your routine testing.

FMPTS

Food Microbiology Proficiency Testing Scheme


Preparation Instruction sheet

Protocol FM2V1M	Sample codes: PT-FM-096, PT-FM-094, PT-FM-092	Description: 2 vials +50g milk powder	
----------------------------	--	---	--


Step 1: For each vial (A+B) weight out 25g of milk powder.



Step 2: Prepare 2*225ml of suitable microbiological diluent.



Step 3: For each vial take 10ml of the 225ml diluent and add it to the vial after aseptically removing cap and rubber stopper.



FMPTS

Food Microbiology Proficiency Testing Scheme

Preparation Instruction sheet


Step 4: Replace the stopper and shake to dissolve and leave the test material to resuscitate (45 min resuscitation).



Step 5: Reconstitute the 25g matrix with the remaining amount of diluent prepared in step 2



Step 6: After resuscitation add the vial contents to the matrix preparation (step5) and do back-washing 2-3 times.




Step 7: Proceed with your routine testing.


FMPTS

Food Microbiology Proficiency Testing Scheme


Preparation Instruction sheet

Protocol FM2V	Sample codes: PT-FM-093, PT-FM-091	Description: 2 vials	
--------------------------	---------------------------------------	-------------------------	---

Step 1: Prepare 2*225ml of suitable microbiological diluent.



Step 3: For each vial take 10ml of the 225ml diluent and add it to the vial after aseptically removing cap and rubber stopper.



FMPTS

Food Microbiology Proficiency Testing Scheme

Preparation Instruction sheet

Step 4: Replace the stopper and shake to dissolve and leave the test material to resuscitate (45 min resuscitation).



Step 7: Proceed with your routine testing.