



# WMPTS-25

## Water Microbiology Proficiency Testing Scheme

# WMPTS

## Water Microbiology Proficiency Testing Scheme Preparation Instruction sheet

Sample ID	Protocol	Page
<u>General Information</u>		3
PT-WM-01	<a href="#">WM-1V</a>	4-6
PT-WM-02		
PT-WM-03		
PT-WM-04		
PT-WM-05		
PT-WM-06		
PT-WM-07		
PT-WM-08		

# **WMPTS**

## **Water Microbiology Proficiency Testing Scheme**

### **Preparation Instruction sheet**

## **GENERAL INFORMATION**

### **Receipt and Storage**

- On receipt of the test material, record the date and store refrigerated at 2-8 °C until ready to test.
- The test material should be analyzed in accordance with the deadlines shown on the application form.
- The test materials represent 'real' water samples, which may or may not contain the target organism(s), at a range of inoculum levels. Test materials may also contain background flora.
- Please consider all dilution factors when calculating the level of the target organisms in the original sample. Where relevant guidance is given in the instructions.
- The test materials should be treated in the same way as routine samples and should therefore be tested using your routine laboratory equipment, reagents and methods.
- Where the instruction mentions a **resuscitate** step, the lyophilized test material should be left at room temperature for between 45 minutes before proceeding with the test.
- Do not sub-divide samples prior to dilution unless stated in the instructions.
- Please treat the sample as contaminated material and use aseptic techniques throughout the test procedure.
- Please ensure samples are mixed well at all stages before proceeding with testing.
- Reconstitute the test material according to the appropriate instructions for each sample type.

### **Precautions**

- Microbiological test materials contain viable micro-organisms and are supplied on the understanding that the purchaser has suitably competent and qualified personnel to handle them safely. Test materials must only be opened in a laboratory by qualified personnel.
- Refer to the Safety Data Sheet for information on the safe handling and disposal of the test materials.

### **Reporting Results**

- All results should be submitted using NLPTS reporting form.

# **WMPTS**

## **Water Microbiology Proficiency Testing Scheme** **Preparation Instruction sheet**


If you need any help at all, please do not hesitate to contact our support team using the details below.

**Tel:** +962 6 5506539

**Email:** [pt-coordinator @naratechlab.com](mailto:pt-coordinator@naratechlab.com)

# WMPTS




## Water Microbiology Proficiency Testing Scheme Preparation Instruction sheet

<b>Protocol WM1V</b>	<b>Sample code:</b>		<b>Description: 1 vial</b>
	PT-WM-01	PT-WM-02	
	PT-WM-03	PT-WM-04	
	PT-WM-05	PT-WM-06	
	PT-WM-07	PT-WM-08	

**Step 1:** prepare 1L sterile water to use as diluent.



**Step 2:** Aseptically remove cap and stopper from the vial and reconstitute the test material by adding 10 ml of the diluent.



# **WMPTS**

## **Water Microbiology Proficiency Testing Scheme**

### **Preparation Instruction sheet**

**Step 3:** replace the vial stopper and shake gently to dissolve.



**Step 4:** leave to stand for 45 min (resuscitation).



**Step 5:** Add this concentrate to the remaining diluent.



# **WMPTS**

## **Water Microbiology Proficiency Testing Scheme**

### **Preparation Instruction sheet**

**Step 6:** repeat this procedure (adding 10 ml of the diluent to the vial and then re add the vial content to the 1L diluent) 2-3 times to ensure all the freeze -dried test material is recovered from the vial.



**Step 7:** invert sample 5 times. (This represents your sample)



**Step 8:** test for the target organism using your routine laboratory method.